

AIV PowerMATE® 20 amp

Technical Specifications

Part Number

PS12595

Dimensions

height - 5.45"

width - 5.28"

depth - 3.27"

weight - 4 lbs

Receptacles

type - Hospital Grade

style - US NEMA 5-20R Hospital Grade

configuration - 2 pole, 3 wire grounding

color - Gray

Circuit Breaker

rating - 20 amps, 125 AC volts

System Rating

15 amps, 120 volts

Power Cord

Hospital Grade

length - 15'

gauge - 12 AWG

plug - US NEMA 5-20P

hospital grade molded plug

color - Black

Pole diameter minimum and maximum

min - 0.875"

max - 1.50"

UL recognized component to UL 1363A



Contact us
aivsales@aiv-inc.com
800.990.2911



FM 86787