

# AIV PowerMATE® 15A - CM

## Technical Specifications

### Part Number

PS13875

### Dimensions

Height - 6.87"

Width - 5.34"

Depth - 3.78"

Weight - 3.25 lbs

### Receptacles

Type - Hubbell Hospital Grade

Style - US NEMA 5-15R Hospital Grade

Configuration - 2 Pole, 3 Wire Grounding

Installed - Ground Pin Up

Color - Gray

### Dual Circuit Breakers

Rating - 15 amps, 125 AC volts

### System Rating

12 amps, 120 volts

### Current Monitoring LED Display

Displays current usage (Green LED)

Power ON indicator (Blue LED)

### Current Bar

Located under current monitoring label

Can be used with a clamp meter to test current usage

### Locking Bolt with Nut and (6) Locking Outlet Covers included

### Power Cord

Over-Molded Hospital Grade

Length - 15'

Gauge - 12 AWG

Plug - US NEMA 5-15P

hospital grade molded plug

Color - Gray

### Pole diameter minimum and maximum

Min - 0.875"

Max - 1.50"

### UL recognized component to UL 1363A

UL recognized 60601-1

CAN/CSA C22.2 60601-1



## Contact us

[aivsales@aiv-inc.com](mailto:aivsales@aiv-inc.com)

800.990.2911



FM 86787